


Lifeline™ Emergency Cart — 35.01

Features include:

- Lightweight polymer and 1 1/4" (32mm) tubular steel construction.
- Ergonomically-placed handle insures maximum control.
- Convertible caster system provides positive control in transit and swivel mode for ease of positioning at code site.
- Acrylic storage areas are designed for multi-person access.
- Recessed top storage area provides quick access to first-line meds.
- Flip-open bins provide secure storage for gloves, masks, etc.
- Smooth-sliding, ball bearing drawer slides provide full extension for complete access to drawer contents.
- Lower storage area for bulk supplies.
- Quick-release latch secures contents until needed.
- Polyethylene body, drawers, steering handle.
- 5" (127mm) swivel casters with resilient rubber tires, one with brake.

Features with Standard Models

Dimensions with LEC steering arm down —
35"W.x35 3/4"H.x22 1/2"D. (889x908x572mm).



*Maximum of 18" (457mm) of any drawer combination, allows 8 1/2" (216mm) bulk storage with 4 1/4" (108mm) opening.

Features	Cat. No. LEC31	Cat. No. LEC41	Cat. No. LEC42	Cat. No. LEC51	Cat. No. LEC52	Cat. No. LEC53	Cat. No. LEC63	Cat. No. LEC75	Cat. No. LEC83
Backboard	X	X	X	X	X	X	X	X	X
Oxygen Tank Holder	X	X	X	X	X	X	X	X	X
I.V. Pole Holder	X	X	X	X	X	X	X	X	X
I.V. Pole			X		X	X	X	X	X
Small Drawer (3") (76mm)	1	1	1	2	2	2	3	4	6
Large Drawer (6") (152mm)	1	1	1	1	1	1	1	1	
Top Cavity Tray						X	X	X	X
Plastic I.V. Bin (Set of 3)		X	X	X	X	X	X	X	X
Defibrillator Swing Arm (Adjustable)			†		†	†	†	†	†
Waste Bag Holder	X								
Plastic Security Seals (Pkg. of 100)	X	X	X	X	X	X	X	X	X
Gel Bin						X	X	X	X
Vacuum Pump/Ambu Bracket						X	X	X	X
Wire Shelf		X	X						
List Price Each	1,630.50	1,772.50	2,173.00	1,852.50	2,253.00	2,397.50	2,513.50	2,768.50	2,882.50



Exchange Tray
LEC9460 (with extra divider) or
LEC9461 (with extra divider)

Accessories

Description	Cat. No.	List Price Each
Adjustable Defibrillator	LEC9003*	343.50
Swing Arm Assembly		
I.V. Pole	LEC9100	93.50
Backboard	LEC9150	136.50
Plastic I.V. Bins (Set of 3)	LEC9200M	181.00
Plastic I.V. Bin (1)	LEC9205M	62.00
Gel Bin	LEC9225	47.50
6" (152mm) Drawer	LEC9300	145.00
3" (76mm) Drawer	LEC9350	129.50
Top Cavity Tray approx.	LEC9460	126.00
15"x15" (381x381mm) with 4 Dividers, 2 ea. front to rear rails, Sheet of Blank Labels, 5 Security Bags		
3" (76mm) Drawer Tray approx.	LEC9461	129.50
20"x15" (508x381mm) with 4 Dividers, 2 ea. front to rear rails, Sheet of Blank Labels, 5 Security Bags		
6" (152mm) Drawer Tray approx. 19"x15" (483x381mm) with Dividers (2 long, 6 short); 5 Security Bags		
6" (152mm) Drawer Divider Kit; 2 Long, 6 Short, 1 Sheet of Blank Labels	LEC9466	156.50

*Defibrillator shelf adjusts 14 1/2"x 9 3/4" (368x248mm) to 14 1/2"x 13 1/8" (368x333mm).

Description	Cat. No.	List Price Each
Clear Lid for LEC9465	LEC9467	55.00
Short Divider for LEC9465 (3 pack)	LEC9468	15.00
Long Divider for LEC9465 (2 pack)	LEC9469	34.00
Security Bags (Seals Over LEC9460/1 for Passive Security, 20/pack)	LEC9470	36.00
Ampule Insert and 12 Dividers for Top Cavity or 3" (76mm) Drawer Trays	LEC9471	89.50
Dividers for LEC9471, 10 pack	LEC9472	41.00
Short Dividers for LEC9460 and LEC9461 10/pack	LEC9473	41.00
Short Dividers for LEC9300 (6" [152mm] Drawer) 5/pack	LEC9474	47.50
Hard Lid for LEC9460	LEC9476	50.50
Long Divider Rails for LEC9460 and LEC9461 (2 pack)	LEC9478	29.50
Long Divider Rails for LEC9300 (6" [152mm] Drawer) (2 pack)	LEC9479	39.00
Vacuum Pump/Ambu Bracket	LEC9600A	82.00
Wire Shelf	LEC9700	49.50
Sharps Container Bracket	LEC9800	42.00
Plastic Security Seals (Pack of 100)	LEC9900A	49.50
Wastebag Holder	LEC-WBH	88.50
†Hospital Grade 6-Outlet Strip 15'	MFLEC595	228.00

†Includes cord wrap.



LEC9205M



LEC9600A

Lifeline Carts are built to exact customer specification, therefore, they are non-cancelable, non-returnable.